

Anti-Hibadh Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31937

Item Number: CM0024245



Introduction

Rabbit polyclonal antibody against human HIBADH, suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Mitochondrial enzyme involved in branched-chain amino acid catabolism.

[Learn More](#)

Parameter	Value
Product Name	HIBADH Rabbit pAb
Remarks/Alias	NS5ATP1; HIBADH
Species	Human
Gene ID (Human)	11112
Gene ID	11112
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-336 of human HIBADH (NP_689953.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31937
Protein Weight	35kDa
Shipping	Ice bag