


Anti-Hibch Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Nvy1

Item Number: CM0020742



Introduction

Rabbit polyclonal antibody against human HIBCH (3-hydroxyisobutyryl-CoA hydrolase), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. 43kDa mitochondrial enzyme in valine catabolism.

[Learn More](#)

Parameter	Value
Product Name	HIBCH Rabbit pAb
Remarks/Alias	HIBYLCOAH; HIBCH
Species	Human
Gene ID (Human)	26275
Gene ID	26275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-330 of human HIBCH (NP_055177.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6NVY1
Protein Weight	43kDa
Shipping	Ice bag