

Anti-Hif1A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q16665

Item Number: CM0011166



Introduction

Rabbit monoclonal antibody against human HIF1 α , validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Suitable for hypoxia research, cancer biology, and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	HIF1 α Rabbit mAb
Remarks/Alias	HIF1; MOP1; PASD8; HIF-1A; bHLHe78; HIF-1alpha; HIF1-ALPHA; HIF-1-alpha
Species	Human
GeneID (Human)	3091
GeneID	3091
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16665
Protein Weight	83 kDa/93 kDa/96 kDa
Shipping	Ice bag