

Anti-Hif1A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16665

Item Number: CM0017740



Introduction

High-quality rabbit polyclonal antibody against HIF1 α (human, mouse, rat) for WB, IHC-P, IF/ICC, ELISA. Targets hypoxia-inducible factor 1-alpha, a master regulator of hypoxia response. Available from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	HIF1 α Rabbit pAb
Remarks/Alias	HIF1; MOP1; PASD8; HIF-1A; bHLHe78; HIF-1alpha; HIF1-ALPHA; HIF-1-alpha; HIF1 α
Species	Human
GeneID (human)	3091
GeneID	3091
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of HIF1 α (NP_001521.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16665
Protein Weight	83 kDa/93 kDa/96 kDa
Shipping	Ice bag