


Anti-Hif1An/Fih-1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nwt6

Item Number: CM0007746



Introduction

Rabbit monoclonal antibody targeting human HIF1AN/FIH-1 (SWISS Q9NWT6). Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Suitable for hypoxia and oxygen sensing research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HIF1AN/FIH-1 Rabbit mAb
Remarks/Alias	FIH1
Species	Human
GeneID (Human)	55662
GeneID	55662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human HIF1AN/FIH1 (NP_060372.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWT6
Protein Weight	41kDa
Shipping	Ice bag