


Anti-Hif1An/Fih1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nwt6

Item Number: CM0024166



Introduction

Rabbit polyclonal antibody targeting HIF1AN/FIH1, validated for Western blot and ELISA. Reacts with human and mouse. Suitable for investigating hypoxia signaling and oxygen sensing in research and assay development.

[Learn More](#)

Parameter	Value
Product Name	HIF1AN/FIH1 Rabbit pAb
Remarks/Alias	FIH1; HIF1AN/FIH1
Species	Human
Gene ID (Human)	55662
Gene ID	55662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human HIF1AN/FIH1 (NP_060372.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NWT6
Protein Weight	41kDa
Shipping	Ice bag