


Anti-Hif1An/Fih1 Rabbit Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Nwt6

Item Number: CM0030414



Introduction

Rabbit polyclonal antibody against HIF1AN (FIH1), validated for knockout studies. Suitable for Western blot, IHC, IF/ICC, IP, ELISA in human, mouse, rat. Targets hypoxia-inducible factor asparagine hydroxylase.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HIF1AN/FIH1 Rabbit pAb
Remarks/Alias	FIH1; H1
Species	Human
GeneID (Human)	55662
GeneID	55662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human HIF1AN/FIH1 (NP_060372.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWT6
Protein Weight	41kDa
Shipping	Ice bag