

Anti-Hif1B/Arnt Polyclonal Antibody For Wb, Elisa - P27540

Item Number: CM0015920



Introduction

Rabbit polyclonal antibody against human HIF1 β /ARNT, suitable for Western blot and ELISA applications. Reacts with human and mouse. Ideal for studying hypoxia and xenobiotic signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	HIF1 β /ARNT Rabbit pAb
Remarks/Alias	HIF1B; TANGO; bHLHe2; HIF1BETA; HIF-1beta; HIF1-beta; HIF-1-beta; HIF1 β /ARNT
Species	Human
GeneID (Human)	405
GeneID	405
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 420-485 of human HIF1 β /ARNT (NP_001659.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27540
Protein Weight	87kDa
Shipping	Ice bag