

Anti-Hif1B/Arnt Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P27540

Item Number: CM0003327



Introduction

Rabbit monoclonal antibody specifically detecting human HIF1 β /ARNT (P27540), validated for WB, IHC-P, and ELISA. Ideal for studying hypoxia signaling and dioxin response pathways.

[Learn More](#)

Parameter	Value
Product Name	HIF1 β /ARNT Rabbit mAb
Remarks/Alias	HIF1B; TANGO; bHLHe2; HIF1BETA; HIF-1beta; HIF1-beta; HIF-1-beta; HIF1 β /ARNT
Species	Human
GeneID	405
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 499-789 of human HIF-1 beta (P27540).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P27540
Protein Weight	87kDa
Shipping	Ice bag