


Anti-Hif3A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2N7

Item Number: CM0023804



Introduction

Rabbit polyclonal antibody against HIF3A (Hypoxia-inducible factor 3-alpha) for Western blot, IF/ICC, ELISA. Validated in human, mouse, rat. Ideal for hypoxia and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	HIF3A Rabbit pAb
Remarks/Alias	IPAS; MOP7; PASD7; HIF-3A; bHLHe17; HIF3-alpha-1; HIF3A
Species	Human
Gene ID (Human)	64344
Gene ID	64344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human HIF3A (NP_690008.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2N7
Protein Weight	72kDa
Shipping	Ice bag