


Anti-Higd1A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y241

Item Number: CM0016400



Introduction

Rabbit polyclonal antibody against human HIGD1A (HIG1, RCF1a), validated for WB, IHC-P, IF/ICC, ELISA. Targets mitochondrial complex IV subunit involved in respiratory chain. Suitable for hypoxia and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	HIGD1A Rabbit pAb
Remarks/Alias	HIG1; RCF1a; HIGD1A
Species	Human
GeneID (Human)	25994
GeneID	25994
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-93 of human HIGD1A (NP_054775.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y241
Protein Weight	10kDa
Shipping	Ice bag