

Anti-Hilpda Polyclonal Antibody For Ihc-P, Elisa - Q9Y5L2

Item Number: CM0022479

Introduction

Rabbit polyclonal antibody against human HILPDA (HIG2) for immunohistochemistry (IHC-P) and ELISA applications. Cross-reacts with human, mouse, and rat. Suitable for hypoxia and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HILPDA Rabbit pAb
Remarks/Alias	HIG2; HIG-2; C7orf68; HILPDA
Species	Human
GeneID (Human)	29923
GeneID	29923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human HILPDA (NP_037464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5L2
Protein Weight	7kDa
Shipping	Ice bag