


Anti-Hint1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49773

Item Number: CM0002528



Introduction

HINT1 Rabbit mAb (P49773) is a monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. It detects HINT1 in human, mouse, and rat samples. Ideal for studying nucleotide metabolism and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	HINT1 Rabbit mAb
Remarks/Alias	HINT, NMAN, PKCI-1, PRKCNH1, HINT1
Species	Human
GeneID (human)	3094
GeneID	3094
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HINT1 (P49773)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49773
Protein Weight	14kDa
Shipping	Ice bag