


Anti-Hint1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49773

Item Number: CM0018774



Introduction

Rabbit polyclonal antibody against human HINT1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Recognizes a 14kDa protein involved in nucleotide metabolism and transcriptional regulation. Suitable for research in apoptosis and cancer.

[Learn More](#)

Parameter	Value
Product Name	HINT1 Rabbit pAb
Remarks/Alias	HINT; NMN; PKCI-1; PRKCNH1; HINT1
Species	Human
GeneID (Human)	3094
GeneID	3094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human HINT1 (NP_005331.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49773
Protein Weight	14kDa
Shipping	Ice bag