

Anti-Hip1 Polyclonal Antibody For Wb And Elisa - 000291

Item Number: CM0029490



Introduction

Rabbit polyclonal antibody targeting human HIP1, ideal for WB and ELISA. Detects human and mouse HIP1. Supports research in clathrin-mediated endocytosis, neurodegeneration, and cancer.

[Learn More](#)

Parameter	Specification
Product Name	HIP1 Rabbit pAb
Remarks/Alias	SHON, HIP-I, ILWEQ, SHONbeta, SHONGgamma, HIP1
Species	Human
GeneID	3092
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-690 of human HIP1 (NP_005329.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00291
Protein Weight	116kDa
Shipping	Ice bag