

Anti-Hip1R Polyclonal Antibody For Wb, If/Icc, Elisa - 075146

Item Number: CM0016015



Introduction

Rabbit polyclonal antibody targeting HIP1R (HIP12), validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat. Ideal for studying clathrin-mediated endocytosis and receptor trafficking.

[Learn More](#)

Parameter	Value
Product Name	HIP1R Rabbit pAb
Remarks/Alias	HIP3; HIP12; ILWEQ; HIP1R
Species	Human
GeneID (human)	9026
GeneID	9026
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human HIP1R (NP_003950.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075146
Protein Weight	119kDa
Shipping	Ice bag