

## Anti-Hipk2 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9H2X6

Item Number: CM0004489

### Introduction

Rabbit monoclonal antibody against human HIPK2, validated for WB, IF/ICC, IF-P, and ELISA. Reacts with human and mouse. Target involved in transcription, apoptosis, and DNA damage response. For research use.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | HIPK2 Rabbit mAb                                                                                                      |
| Remarks/Alias    | PRO0593; hHIPk2; HIPK2                                                                                                |
| Species          | Human                                                                                                                 |
| GeneID (Human)   | 28996                                                                                                                 |
| GeneID           | 28996                                                                                                                 |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HIPK2 (Q9H2X6). |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                                          |
| SWISS            | Q9H2X6                                                                                                                |
| Protein Weight   | 101KDa/128KDa/131kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                               |