


Anti-Hippocalcin (Hpca) Rabbit Monoclonal Antibody For Wb, Elisa - P84074

Item Number: CM0002286



Introduction

High-quality rabbit monoclonal antibody against human Hippocalcin (HPCA), validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for neuroscience research and calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Hippocalcin (HPCA) Rabbit mAb
Remarks/Alias	BDR2; DYT2; Hippocalcin (HPCA)
Species	Human
Gene ID (Human)	3208
Gene ID	3208
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Hippocalcin (HPCA) (P84074)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84074
Protein Weight	22kDa
Shipping	Ice bag