


Anti-Hira Monoclonal Antibody For Wb, Elisa - P54198

Item Number: CM0011821



Introduction

Rabbit monoclonal antibody against human HIRA, validated for Western blot and ELISA. Protein HIRA is a chromatin regulator involved in replication-independent chromatin assembly and cellular senescence. Ideal for research in cell cycle and chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	HIRA Rabbit mAb
Remarks/Alias	TUP1, DGCR1, TUPLE1, HIRA
Species	Human
GeneID (Human)	7290
GeneID	7290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-675 of human HIRA (P54198)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P54198
Protein Weight	112kDa
Shipping	Ice bag