


Anti-Hira Polyclonal Antibody For Wb And Elisa - P54198

Item Number: CM0017907



Introduction

Rabbit polyclonal antibody against human HIRA (TUP1/DGCR1), a chromatin assembly factor and transcriptional repressor. Validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for epigenetics, cell cycle, and senescence research.

[Learn More](#)

Parameter	Value
Product Name	HIRA Rabbit pAb
Remarks/Alias	TUP1, DGCR1, TUPLE1, HIRA
Species	Human
GeneID (human)	7290
GeneID	7290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human HIRA (NP_003316.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54198
Protein Weight	112kDa
Shipping	Ice bag