

Anti-Hirip3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bw71

Item Number: CM0028934

Introduction

High-quality rabbit polyclonal antibody against human HIRIP3, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Histone chaperone facilitating H2A-H2B deposition onto chromatin.

[Learn More](#)

Parameter	Value
Product Name	HIRIP3 Rabbit pAb
Remarks/Alias	HIRIP3
Species	Human
Gene ID	8479
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 417-556 of human HIRIP3 (NP_003600.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BW71
Protein Weight	62kDa
Shipping	Ice bag