


## Anti-Hist1H3D Polyclonal Antibody For Wb, Elisa - Q16695 / P68431

Item Number: CM0023055



### Introduction

Rabbit polyclonal antibody targeting HIST1H3D (H3.4/H3/g), validated for Western blot and ELISA. Predicted to react with human, mouse, rat, and other species. Ideal for chromatin biology and histone modification studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | HIST1H3D Rabbit pAb                                                                                                  |
| Remarks/Alias    | H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; HIST1H3D                                                       |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 8290/8350                                                                                                            |
| Gene ID          | 8290/8350                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human HIST1H3D (NP_003520.1). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                                      |
| SWISS            | Q16695/P68431                                                                                                        |
| Protein Weight   | 15KDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |