


Anti-Histamine H2 Receptor Polyclonal Antibody For Wb And Elisa - P25021

Item Number: CM0016255



Introduction

Rabbit polyclonal antibody targeting the human histamine H2 receptor (HRH2). Validated for Western blot and ELISA applications. Detects human, mouse, and rat samples. Suitable for gastric acid secretion and GPCR signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HRH2 Rabbit pAb
Remarks/Alias	H2R; HH2R; HRH2
Species	Human
Gene ID (Human)	3274
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human HRH2 (NP_001124527.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25021
Protein Weight	40kDa
Shipping	Ice bag