

Anti-Histamine H3 Receptor (Hrh3) Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Y5N1

Item Number: CM0004030



Introduction

Rabbit monoclonal antibody against human Histamine H3 receptor (HRH3, GPCR97). Validated for WB, IP, ELISA. Reacts with human, mouse, rat. Suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	HRH3 Rabbit mAb
Remarks/Alias	HH3R; GPCR97; HRH3
Species	Human
GeneID (Human)	11255
GeneID	11255
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HRH3 (Q9Y5N1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5N1
Protein Weight	49kDa
Shipping	Ice bag