

Anti-Histamine H3 Receptor (Hrh3) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5N1

Item Number: CM0019699



Introduction

High-quality rabbit polyclonal antibody against human histamine H3 receptor (HRH3), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and GPCR research.

[Learn More](#)

Parameter	Value
Product Name	HRH3 Rabbit pAb
Remarks/Alias	HH3R, GPCR97, HRH3
Species	Human
GeneID	11255
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human HRH3 (NP_009163.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5N1
Protein Weight	49kDa
Shipping	Ice bag