


Anti-Histamine H4 Receptor (Hrh4) Polyclonal Antibody For Wb And Elisa - Q9H3N8

Item Number: CM0016038



Introduction

Rabbit polyclonal antibody targeting human Histamine H4 receptor (HRH4). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Designed for studying immune cell chemotaxis and inflammatory signaling research.

[Learn More](#)

Parameter	Value
Product Name	HRH4 Rabbit pAb
Remarks/Alias	H4; H4R; BG26; HH4R; AXOR35; GPRv53; GPCR105; HRH4
Species	Human
Gene ID (Human)	59340
Gene ID	59340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 194-304 of human HRH4 (NP_067637.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3N8
Protein Weight	44kDa
Shipping	Ice bag