

Anti-Histidine Ammonia-Lyase (Hal) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42357

Item Number: CM0017754

Introduction

Rabbit polyclonal antibody targeting Histidine ammonia-lyase (HAL). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for histidine metabolism and disease variant studies.

[Learn More](#)

Parameter	Value
Product Name	HAL Rabbit pAb
Remarks/Alias	HIS; HSTD; HAL
Species	Human
GeneID (Human)	3034
GeneID	3034
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HAL (NP_002099.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42357
Protein Weight	73kDa
Shipping	Ice bag