

Anti-Histidine-Rich Glycoprotein (Hrg) Polyclonal Antibody For Wb, Elisa - P04196

Item Number: CM0026174

Introduction

Rabbit polyclonal antibody targeting human Histidine-rich glycoprotein (HRG). Validated in WB and ELISA, cross-reacts with mouse and rat. Ideal for studying coagulation, angiogenesis, and immune clearance mechanisms.

[Learn More](#)

Parameter	Value
Product Name	HRG Rabbit pAb
Remarks/Alias	HPRG; HRGP; THPH11; HRG
Species	Human
GeneID (Human)	3273
GeneID	3273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human HRG (NP_000403.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04196
Protein Weight	60kDa
Shipping	Ice bag