


Anti-Histone H1.0 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07305

Item Number: CM0009965



Introduction

Rabbit monoclonal antibody targeting human Histone H1.0 (P07305), validated for Western blot, IHC-paraffin, and ELISA; cross-reacts with mouse and other species; 21kDa; ideal for chromatin biology and differentiation research.

[Learn More](#)

Parameter	Specification
Product Name	Histone H1.0 Rabbit mAb
Remarks/Alias	H10; H1.0; H1F0; H1FV; Histone H1.0
Species	Human
GeneID (Human)	H1
GeneID	3005
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H1.0 (P07305).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Other (Wide Range Predicted)
SWISS	P07305
Protein Weight	21kDa
Shipping	Ice bag