


Anti-Histone H1.2 (Hist1H1C) Rabbit Polyclonal Antibody - P16403

Item Number: CM0028253



Introduction

Rabbit polyclonal antibody targeting human Histone H1.2 (HIST1H1C). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for chromatin condensation and epigenetic studies.

[Learn More](#)

Parameter	Specification
Product Name	HIST1H1C Rabbit pAb
Remarks/Alias	H1C; H1.2; H1F2; H1s-1; HIST1H1C
Species	Human
GeneID (Human)	3006
GeneID	3006
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16403
Protein Weight	21kDa
Shipping	Ice bag