


Anti-Histone H1.2 Monoclonal Antibody For Wb, If/Icc, Elisa - P16403

Item Number: CM0002359



Introduction

Rabbit monoclonal antibody against human Histone H1.2 (P16403). Validated for WB, IF/ICC, ELISA. Predicted wide cross-reactivity. Suitable for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Histone H1.2 Rabbit mAb
Remarks/Alias	H1C; H1.2; H1F2; H1s-1; HIST1H1C; Histone H1.2
Species	Human
GeneID (Human)	3006
GeneID	3006
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H1.2 (P16403))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Other (Wide Range Predicted)
SWISS	P16403
Protein Weight	21kDa
Shipping	Ice bag