


Anti-Histone H1.4 Rabbit Polyclonal Antibody For Wb And Elisa - P10412

Item Number: CM0019349



Introduction

Rabbit polyclonal antibody targeting human Histone H1.4, suitable for western blot and ELISA with cross-reactivity to human, mouse, and rat. Ideal for chromatin biology and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Histone H1.4 Rabbit pAb
Remarks/Alias	H1E; H1.4; H1F4; RMNS; H1s-4; HIST1H1E; dj221C16.5; Histone H1.4
Species	Human
GeneID (Human)	3008
GeneID	3008
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H1.4 (NP_005312.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10412
Protein Weight	22kDa
Shipping	Ice bag