

Anti-Histone H1T Polyclonal Antibody For Wb, Ihc-P, Elisa - P22492

Item Number: CM0019971

Introduction

Rabbit polyclonal antibody targeting human Histone H1t, a testis-specific linker histone. Confirmed for WB, IHC-P, ELISA with reactivity across human, mouse, rat. Suitable for chromatin biology and spermatogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Histone H1t Rabbit pAb
Remarks/Alias	H1t; H1.6; H1FT; HIST1H1T; dj221C16.2
Species	Human
GeneID (Human)	H1
GeneID	3010
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-207 of human Histone H1t (NP_005314.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22492
Protein Weight	22kDa
Shipping	Ice bag