


## Anti-Histone H2A Monoclonal Antibody For Wb, Ihc-P, Elisa - P04908/P0C0S8

Item Number: CM0008546



### Introduction

Rabbit monoclonal antibody targeting human Histone H2A (P04908/P0C0S8). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat, and predicted broad species. Ideal for chromatin and epigenetics research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Histone H2A Rabbit mAb                                                                                                    |
| Remarks/Alias    | H2A.1; H2A.2; H2A/a; H2AC4; H2AFA; HIST1H2AE; Histone H2A                                                                 |
| Species          | Human                                                                                                                     |
| GeneID (human)   | 3012                                                                                                                      |
| GeneID           | 3012/8329                                                                                                                 |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H2A (P04908). |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | WB, IHC-P, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                                           |
| SWISS            | P04908/P0C0S8                                                                                                             |
| Protein Weight   | 14 kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                   |