


Anti-Histone H2Afv Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q71Ui9

Item Number: CM0019487



Introduction

CamelBio offers anti-Histone H2AFV (H2A.Z-2) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Histone H2AFV Rabbit pAb
Remarks/Alias	H2AV; H2AFV; H2A.Z-2; Histone H2AFV
Species	Human
GeneID (Human)	94239
GeneID	94239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human Histone H2AFV (NP_036544.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q71UI9
Protein Weight	13kDa
Shipping	Ice bag