


Anti-Histone H2Ax Monoclonal Antibody For Ihc-P, Chip, Elisa - P16104

Item Number: CM0002729



Introduction

High-specificity rabbit monoclonal antibody targeting Histone H2AX. Suitable for IHC-P, ChIP, and ELISA. Cross-reacts with human, mouse, rat, and more. Ideal for DNA damage and chromatin research.

[Learn More](#)

Parameter	Specification
Product Name	Histone H2AX Rabbit mAb
Remarks/Alias	H2A.X; H2A/X; H2AFX; Histone H2AX
Species	Human
GeneID (human)	H2A
GeneID	3014
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H2AX (P16104).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P16104
Protein Weight	15kDa
Shipping	Ice bag