


Anti-Histone H2Ax Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16104

Item Number: CM0018439



Introduction

Rabbit polyclonal antibody against Histone H2AX (Human, Mouse, Rat) validated for WB, IHC-P, IF/ICC and ELISA. Suitable for DNA damage response and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	Histone H2AX Rabbit pAb
Remarks/Alias	H2A.X; H2A/X; H2AFX; Histone H2AX
Species	Human
Gene ID (Human)	H2A
Gene ID	3014
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P16104
Protein Weight	15kDa
Shipping	Ice bag