


Anti-Histone H2A.z Monoclonal Antibody For Wb, Ip, Elisa, Chip - P0C0S5

Item Number: CM0004786



Introduction

Rabbit monoclonal antibody against human Histone H2A.Z, validated by knockout, for WB, IP, ELISA, ChIP. Cross-reacts with human, mouse, rat, and other species. Ideal for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Histone H2A.Z Rabbit mAb
Remarks/Alias	H2AZ; H2A.z; H2A/z; H2AFZ; H2A.Z-1; .Z
Species	Human
GeneID (Human)	3015
GeneID	3015
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H2A.Z (P0C0S5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P0C0S5
Protein Weight	14kDa
Shipping	Ice bag