


Anti-Histone H2B Polyclonal Antibody For Wb, Ihc-P, Elisa - P62807 / Q16778

Item Number: CM0023630



Introduction

Rabbit polyclonal antibody against human Histone H2B (P62807/Q16778) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat, and wide range predicted. Suitable for chromatin and nucleosome research.

[Learn More](#)

Parameter	Value
Product Name	Histone H2B Rabbit pAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; Histone H2B
Species	Human
GeneID (Human)	H2B
GeneID	3017/8349
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Histone H2B (NP_003519.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807/Q16778
Protein Weight	14kDa
Shipping	Ice bag