


Anti-Histone H2B Rabbit Monoclonal Antibody For Wb - P62807 / Q16778

Item Number: CM0005361



Introduction

Rabbit monoclonal antibody against human Histone H2B, validated for WB, IHC-P, ChIP, ELISA. Cross-reacts with human, mouse, rat. Ideal for chromatin and epigenetics research applications.

[Learn More](#)

Parameter	Value
Product Name	Histone H2B Rabbit mAb
Remarks/Alias	H2B; H2BE; H2BQ; GL105; H2B.1; H2BFQ; H2BGL105; H2B-GL105; HIST2H2BE; Formyl-Histone H2B-K108
Species	Human
GeneID (Human)	3017/8349
GeneID	3017/8349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human Histone H2B (NP_003519.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62807/Q16778
Protein Weight	14 kDa
Shipping	Ice bag