


Anti-Histone H2B Rabbit Polyclonal Antibody For Wb, Elisa - P62807 / Q16778

Item Number: CM0020713



Introduction

Anti-Histone H2B rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Recognizes human, mouse, rat, and other species with predicted wide cross-reactivity. Ideal for chromatin, nucleosome, and PTM research.

[Learn More](#)

Parameter	Details
Product Name	Histone H2B Rabbit pAb
Remarks/Alias	H2B, H2BE, H2BQ, GL105, H2B.1, H2BFQ, H2BGL105, H2B-GL105, HIST2H2BE, Histone H2B
Species	Human
GeneID (human)	H2B
GeneID	3017, 8349
Immunogen	Synthetic peptide (aa 1-100 of human Histone H2B, NP_003519.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62807 / Q16778
Protein Weight	14 kDa
Shipping	Ice bag