


Anti-Histone H2B (Testis Specific) Monoclonal Antibody For Wb, Ip, Ihc - Q96A08

Item Number: CM0006770



Introduction

Rabbit monoclonal antibody against human testis-specific Histone H2B type 1-A (H2BC1/TSH2B), validated for WB, IP, IF, IHC, ELISA. Cross-reacts with human, mouse, rat, and wide range predicted. Ideal for chromatin and spermatogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	Histone H2B (testis specific) Rabbit mAb
Remarks/Alias	STBP; TH2B; H2BFU; TSH2B; hTSH2B; TSH2B.1; HIST1H2BA; bA317E16.3
Species	Human
GeneID (human)	255626
GeneID	255626
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Histone H2B (testis specific) (NP_733759.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q96A08
Protein Weight	14kDa
Shipping	Ice bag