


Anti-Histone H3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa, Chip - Q16695 / P68431

Item Number: CM0023484



Introduction

High-quality rabbit polyclonal antibody against Histone H3 (H3.1t/H3.4), validated for WB, IHC-P, IP, ELISA, ChIP. Cross-reacts with human, mouse, rat and other species. Suitable for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Histone H3 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Histone H3
Species	Human
Gene ID (Human)	H3
Gene ID	8290/8350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-136 of human HIST3H3 (NP_003484.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag