


Anti-Histone H3 Rabbit Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q16695 / P68431

Item Number: CM0024391



Introduction

Rabbit polyclonal antibody targeting human Histone H3.1t (H3.4, H3/g). Suitable for WB, IP, IF-P, IHC-P, ELISA. Cross-reacts with mouse, rat, and other species. For chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Histone H3 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Histone H3
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide; a synthetic peptide corresponding to a sequence within amino acids 36-136 of human HIST3H3 (NP_003520.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag