

## Anti-Histone H3.3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa, Chip, Chip-Seq - P84243

Item Number: CM0004860

### Introduction

High-quality anti-Histone H3.3 rabbit monoclonal antibody validated for WB, IHC-P, ELISA, ChIP, and ChIP-seq. Suitable for human, mouse, rat, and other species. Ideal for chromatin and epigenetic research.

[Learn More](#)

| Parameter        | Value                                                                                |
|------------------|--------------------------------------------------------------------------------------|
| Product Name     | Histone H3.3 Rabbit mAb                                                              |
| Remarks/Alias    | H3F3; H3-3B; H3.3A; H3F3A; BRYLIB1; Histone H3.3                                     |
| Species          | Human                                                                                |
| GeneID (Human)   | H3                                                                                   |
| GeneID           | 3020                                                                                 |
| Immunogen        | Synthetic peptide corresponding to amino acids 57-136 of human Histone H3.3 (P84243) |
| Source           | Rabbit                                                                               |
| Category         | Monoclonal Antibodies                                                                |
| Application      | WB, IHC-P, ELISA, ChIP, ChIP-seq                                                     |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                      |
| SWISS            | P84243                                                                               |
| Protein Weight   | 15kDa                                                                                |
| Shipping         | Ice bag                                                                              |