


Anti-Histone H3.3 Rabbit Polyclonal Antibody For Wb, Elisa - Q16695 / P68431

Item Number: CM0018497



Introduction

Rabbit polyclonal antibody against Histone H3.3 (H3F3A/B), validated for WB and ELISA, with wide cross-reactivity predicted across human, mouse, rat, and other species. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Histone H3.3 Rabbit pAb
Remarks/Alias	H3F3; H3-3B; H3.3A; H3F3A; BRYLIB1; Histone H3.3
Species	Human
GeneID (Human)	H3
Gene ID	8290/8350
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human H3F3A (NP_002098.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag