


Anti-Histone H4 Polyclonal Antibody For Wb, Elisa - P62805

Item Number: CM0016504



Introduction

Rabbit polyclonal antibody against Histone H4 (P62805). Validated for WB and ELISA with wide cross-reactivity (human, mouse, rat). Essential tool for chromatin biology and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	Histone H4 Rabbit pAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A; Histone H4
Species	Human
GeneID (Human)	H4
GeneID	8359
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag