


Anti-Histone H4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P62805

Item Number: CM0002868



Introduction

CamelBio offers high-quality Anti-Histone H4 Rabbit mAb validated for WB, IHC-P, ELISA. Detects human, mouse, rat Histone H4 (P62805) with wide cross-reactivity. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Histone H4 Rabbit mAb
Remarks/Alias	H4/p; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; HIST4H4; Histone H4
Species	Human
GeneID (Human)	121504
GeneID	8359
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 24-103 of human Histone H4 (P62805).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag