


Anti-Histone H4 Rabbit Polymab® For Wb, If/Icc, Ihc-P, Chip, Elisa - P62805

Item Number: CM0003549



Introduction

Rabbit recombinant monoclonal antibody (PolymAb®) against human Histone H4. Validated for WB, IF/ICC, IHC-P, ChIP, ELISA. Cross-reacts with mouse, rat. Detects 11 kDa nucleosome core protein for chromatin/epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Histone H4 Rabbit PolymAb®
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A
Species	Human
Gene ID (Human)	8359
Gene ID	8359
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag