


Anti-Hjurp Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ncd3

Item Number: CM0022804



Introduction

HJURP rabbit polyclonal antibody for human samples, validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Essential for centromere and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	HJURP Rabbit pAb
Remarks/Alias	FAKTS; URLC9; hFLEG1; HJURP
Species	Human
GeneID (Human)	55355
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human HJURP (NP_060880.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8NCD3
Protein Weight	84kDa
Shipping	Ice bag